


SURFACE
DISINFECTION



INSTRUMENT
DISINFECTANT

Aldasan® 2000/ Aldosan® 2000-CH

Disinfection of medical devices and other objects
through submersion or wiping



CE0482

Product Information

Please read carefully!

Aldasan® 2000/Aldosan® 2000-CH

Concentrate

For professional use.

INGREDIENTS

100 g of the solution contains:

9.8 g formaldehyde, 9.9 g glutaraldehyde, sodium alkyl ether sulphate, Isotridecanol (ethoxylated)

SCOPE OF APPLICATION

- For the **disinfection** and cleaning in separate steps in case of bacterial (including TBC) or viral infections. Suitable for use on **instruments** and **flexible endoscopes**, as well as inventory and surfaces (medical devices). Only use in the manner indicated in the instructions for use provided by the manufacturer of the medical device.
- For the prevention of nosocomial infections.
- Also suitable for use in ultrasonic baths.

INSTRUCTIONS FOR USE/DOSAGE

Aldasan® 2000/Aldosan® 2000-CH is a liquid concentrate, the desired solution can be reached by adding the requisite amount of water to the concentrate. The water should not exceed room temperature.

Caution: Always first fill the container for the solution with water, and add **Aldasan® 2000/Aldosan® 2000-CH** afterwards. The solution should be prepared immediately before application (in compliance with the directive for hospital hygiene and infection prevention published by the Robert-Koch-Institute, Berlin). Solutions without organic load (i.e., in case of disinfection of pre-cleaned instruments) retain efficacy for at least 14 days.

WARNING AND STORAGE NOTICES

- Keep disinfectant at a cool but frost-free, well ventilated and dry place and out of reach of children.
- Wear protective gloves, protective clothing and eye protection when handling the undiluted concentrate.
- Only use in a well-ventilated environment.
- Ventilate after use on surfaces, or use in a closed container (instrument disinfection).
- Do not use after expiration date.
- Shelf life of solutions in clean conditions: 14 days**.
- Use **Aldasan® 2000/Aldosan® 2000-CH** safely! Always consult label and product information before use.

PRODUCT FEATURES

- bactericidal, yeasticidal, fungicidal, tuberculocidal, mycobactericidal, and virucidal
- Disinfectant that is gentle to the material with comprehensive spectrum of efficacy for use on inventory, textiles, and instruments.
- Excellent cleaning capacity even in dirty conditions. Also suitable for use in the presence of organic substances (e.g., fat, proteins, or serum),
- fresh scent
- Resistant to lime water.
- biodegradable

** The period during which the efficacy of the disinfectant remains preserved, provided the solution was mixed according to the manufacturer's instructions, using the correct dosage, and was left covered and unused during the time frame listed by the manufacturer.

Wiping disinfection acc. to VAH-methods 2015 with mechanical action, dirty conditions

bactericidal, yeasticidal	1.00 % – 15 min
	0.75 % – 30 min
	0.50 % – 1 h
	0.25 % – 4 h

Instrument disinfection acc. to VAH-methods 2015, dirty conditions

bactericidal, yeasticidal	3.00 % – 15 min
	2.00 % – 30 min
	1.00 % – 1 h
mycobactericidal	3.00 % – 30 min
	2.00 % – 1 h

Testing acc. to DVV/RKI Guideline

virucidal activity against enveloped viruses* (incl. HIV/HBV/HCV)	0.50 % – 15 min
adenovirus	1.00 % – 5 min
papovavirus/polyomavirus/SV40	0.50 % – 1 h
poliovirus	2.00 % – 15 min

For prevention of infection acc. to § 18 of the German IfSG (Law for protection against infection)

surface disinfection	4.00 % – 4 h
----------------------	--------------

European Standards (dirty conditions)

EN 13624 fungicidal	1.00 % – 30 min
EN 16615 bactericidal, yeasticidal	1.00 % – 15 min 0.75 % – 30 min 0.50 % – 1 h
EN 14476 limited spectrum of virucidal activity *	1.00 % – 15 min
EN 14476 virucidal activity *	3.00 % – 30 min
EN 17111 virucidal activity *	3.00 % – 30 min

*acc. to recommendations of Robert-Koch-Institute Berlin, published in Federal Health Gazette 3/2017

Caution!

Always place the objects into the disinfecting solution in such a manner that all surfaces are completely submerged. Close container with a lid. After disinfection/cleaning, instruments should be rinsed with water. Endoscopes must be rinsed with sterile water.

CAUTION: The amount of water to be added is calculated by subtracting the requisite amount of disinfectant from the desired solution.

Lysoform · Systematic Hygiene

Concentration	Concentration		
	0.25%	0.50%	0.75%
Dosage System (in 8 litres of water)			
Dosage pump (1 pump = 20 ml/30 ml)	1 pump	2 pumps	4 pumps
Lysoform dosage aid DH 91 – no electronics	One-time calibration for desired concentration		

SPECIAL NOTICES

- a) Product incompatibilities
The addition of household cleaning agents may compromise or neutralize the disinfecting function of **Aldasan® 2000/Aldasan® 2000-CH**.
- b) Countermeasures in case of improper use
Should, through improper use, the product **Aldasan® 2000/Aldasan® 2000-CH** have come into contact with any mucous membranes (e.g., eyes, throat, or stomach), a doctor should be consulted immediately. As an immediate measure, dilution of the preparation is recommended. In case of oral ingestion, drink plenty of water. In case of contact with the eyes, rinse thoroughly with generous amounts of water.
- c) Use **Aldasan® 2000/Aldasan® 2000-CH** safely!
Always consult label and product information before use.
Certified and registered in all relevant listings.

SHIPPING UNIT/CONTAINER

6 x 1 litre bottle
2 x 5 litre canister

Art. no.

6241
6255

Accessories

Dosage pump for 5-litre canister
(1 pump = 20 ml)

3135

Technical aids

Lysoform dosage aid

9174



Lysoform Dr. Hans Rosemann GmbH
Kaiser-Wilhelm-Str. 133, D-12247 Berlin
Telefon +49 30 77992-0
www.lysoform.com, e-mail: kontakt@lysoform.de