

Lysoform Product Catalogue


Products and Accessories




HAND
DISINFECTANTS



SKIN/HAND
PROTECTION/CARE



SURFACE
DISINFECTANTS



INSTRUMENT
DISINFECTANTS



ACCESSORIES



ADDITIONAL
PRODUCTS

**Disinfection and hygiene
since 1900**



Certificate

The Certification Body

MEDCERT Zertifizierungs- und Prüfungsgesellschaft für die Medizin GmbH
Pilotuspool 2 – 20355 Hamburg – Germany

herewith certifies that the company

Lysaform Dr. Hans Rosemann GmbH
Kaiser-Wilhelm-Straße 133
12247 Berlin
Germany

has introduced, applies and maintains a quality management system in the area of activities, products/services and locations listed in the appendix.

The conformity of this quality management system to the requirements of the following standard has been verified by an audit:

EN ISO 13485:2016

This certification is subject to surveillance by MEDCERT.

Effective date: 2022-06-14
Expiry date: 2025-04-30

Report No.: 1208F525F
Procedure No.: QS - 1208
Certificate No.: 1208DE445220614

Hamburg, 2022-06-14

M.H.
MEDCERT Certification Body
Morris Harder

The certificate is only valid when provided entirely with all of its pages.
To verify the validity of this certificate, contact info@medcert.de.

Zertifikat

Die Zertifizierungsstelle

MEDCERT Zertifizierungs- und Prüfungsgesellschaft für die Medizin GmbH
Pilotuspool 2 – 20355 Hamburg – Deutschland

bescheinigt hiermit, dass das Unternehmen

Lysaform Dr. Hans Rosemann GmbH
Kaiser-Wilhelm-Straße 133
12247 Berlin
Deutschland

für den Bereich der in der Anlage genannten Aktivitäten, Produkte/Dienstleistungen und Standorte ein Qualitätsmanagementsystem eingeführt hat, angewendet und aufrechterhält.

Durch ein Audit wurde der Nachweis erbracht, dass das Qualitätsmanagementsystem die Anforderungen der unten genannten Norm erfüllt:

EN ISO 13485:2016

Dies wird von MEDCERT überwacht.

Gültig ab: 2022-06-14
Gültig bis: 2025-04-30

Berichts-Nr.: 1208F525F
Verfahrens-Nr.: QS - 1208
Zertifikats-Nr.: 1208DE445220614

Hamburg, 2022-06-14

M.H.
MEDCERT Zertifizierungsstelle
Morris Harder

EC Certificate of Conformity

The Notified Body

MEDCERT Zertifizierungs- und Prüfungsgesellschaft für die Medizin GmbH
Pilotuspool 2 – 20355 Hamburg – Germany

herewith certifies that the company

Lysaform Dr. Hans Rosemann GmbH
Kaiser-Wilhelm-Straße 133
12247 Berlin
Germany

has introduced and maintains a quality assurance system for the products / product categories listed in the appendix.

The compliance of this quality assurance system with the below mentioned requirements of the Council Directive 93/42/EEC was verified by an audit:

Annex II without section 4

This certification is subject to surveillance by MEDCERT.

This certificate is valid until 30 April 2024

Report No.: 1208F525F
Process No.: QS - 1208
Certificate No.: 1208DE410190801

Hamburg, 01 August 2019

A.S.
MEDCERT Certification Body
Dr. Andreas Schütz

The certificate is only valid when provided entirely with all of its pages.
To verify the validity of this certificate, contact info@medcert.de.

Certification Number: 0482

Rev. 10 / 2019.05.22



EG-Konformitätsbescheinigung

Die Benannte Stelle

MEDCERT Zertifizierungs- und Prüfungsgesellschaft für die Medizin GmbH
Pilotuspool 2 – 20355 Hamburg – Deutschland

bescheinigt hiermit, dass das Unternehmen

Lysaform Dr. Hans Rosemann GmbH
Kaiser-Wilhelm-Straße 133
12247 Berlin
Deutschland

für die Produkte / Produktkategorien gemäß Anlage ein Qualitätssicherungssystem eingeführt hat, angewendet und aufrechterhält.

Durch ein Audit wurde der Nachweis erbracht, dass dieses Qualitätssicherungssystem die unten genannten Anforderungen der Richtlinie 93/42/EWG des Rates erfüllt.

Anhang II ohne Berücksichtigung von Abschnitt 4

Dies wird von MEDCERT überwacht.

Diese Bescheinigung ist gültig bis: 30. April 2024

Berichts-Nr.: 1208F525F
Verfahrens-Nr.: QS - 1208
Bescheinigungs-Nr.: 1208DE410190801

Hamburg, den 01. August 2019

A.S.
MEDCERT Zertifizierungsstelle
Dr. Andreas Schütz

Die Bescheinigung ist nur gültig, wenn sie mit allen Seiten vollständig vorliegt. Bitte wenden Sie sich an info@medcert.de, um die Gültigkeit der Bescheinigung zu überprüfen.

MEDCERT Kennnummer: 0482

Form F10010005 06 / Rev. 10 / 2019.05.22



TABLE OF CONTENTS

Product Overview by Application	4
Product Overview by Main/Active Ingredients and Dosage	5
HyBag®-System	7
Product Overview from A to Z	8
Accessories	16

Product Overview by Application

HANDS

Disinfection

AHD 2000® /AHD 2000® G
Hospisept®
Neoseptin®
ACO®-derm sensitive

Hygienic hand washing

Amosept®

Cleaning

Wasa® 2000
Wasa®-Reinigungsschaum
Wasa®-Soft
ACO®-san sensitive

Care/Protection

Luphenil®
Majola®-H5-Creme
ACO®-care sensitive

SKIN

Desinfektion

AHD 2000®
Hospisept®
Hospisept® (gefärbt)

Cleaning

Wasa® 2000
Wasa®-Reinigungsschaum
Wasa®-Soft
ACO®-san sensitive

Care/protection

Luphenil®
Majola®-H5-Creme
ACO®-care sensitive

MISCELLANEOUS COSMETIC

PRODUCTS

Lysoform® Fuß-Spray

SURFACES/EQUIPEMENT/ INSTRUMENTS

Disinfection with cleaning action

Aldasan® 2000 / Aldosan® 2000 (CH)
Hospisept®-Wipes
Lysoform®
Lysoformin® 3000
Lysoformin® rapid
Lysoformin® rapid wipes
Lysoformin® spezial
Trichlorol®

Disinfection and Cleaning foam

Lysoformin® Plus-Foam

Spray disinfection

Aerodesin® 2000
Fugaten®-Spray

Cleaning

Blanisol®-Pur
Blanisol®-Spezialschaum
LY®-Dry

Rapid disinfection

Aerodesin®
Aerodesin® 2000
Clorina®
Fugaten®
Fugaten®-Spray

Odour eliminator

LY®-Fresh

INSTRUMENTS

Disinfection with cleaning action

Aldasan® 2000 / Aldosan® 2000 (CH)
Almyrol®
Lysoformin® 3000

Cleaning

Blanisol®-Pur

EXCRETIONS

Disinfection

Clorina®
Trichlorol®

LAUNDRY

Disinfection with cleaning action

Trichlorol®

APPLICATION AIDS

Dosage aid DH 1
Dosage aid DH 91
Dosage pumps
Wall-mounted dispenser
Lysoform Soft-Vlies
(container + roll)

Product Overview by Main/Active Ingredients and Dosage

PRODUCT	MAIN INGREDIENTS	*DOSAGE INFORMATION
ACO®-care sensitive	Glyceryl stearate/Isopropyl Palmitate	Ready to use
ACO®-derm sensitive (VAH, IHO)	Ethanol	Ready to use
ACO®-san sensitive	Tensides/refattening agents	Ready to use
Aerodesin® (VAH, IHO)	Propan-1-ol/Ethanol	Ready to use
Aerodesin® 2000 (VAH, IHO)	Propan-1-ol/Ethanol/ glutaraldehyde	Ready to use
AHD 2000®/AHD 2000® G (VAH, RKI, IHO)	Ethanol	Ready to use
Aldasan® 2000 (VAH, RKI, IHO)	Formaldehyde/glutaraldehyde/ anion. Tensides	Surfaces 0.25 % Instruments incl. Mycobact. terrae
Almyrol® (VAH, IHO)	Guanidinderivates/Alkylamin- derivates/Quat	Instruments 0.5 % – 2.0 %
Amosept® (VAH, IHO)	Quat/ampho. Tensides	Ready to use
Blanisol®-Pur	anion. Tensides	0.5 %
Blanisol® Spezialschaum	anion. Tensides/anion. Tensides	Ready to use
Clorina® (VAH, RKI, IHO)	Tosylchloramide sodium	0.5 %
Fugaten®	Ethanol	Ready to use
Fugaten®-Spray (VAH, IHO)	Ethanol	Ready to use
Hospisept® (VAH, RKI, IHO)	Propan-1-ol/Ethanol	Ready to use
Hospisept®-wipes	Propan-1-ol/Ethanol	Ready to use
Luphenil®	Water/Oil emulsion/Beeswax/ Zinkstearate/Caring agents	Ready for use
LY®-Dry		Ready for use
LY®-Fresh		Ready for use

Product Overview by Main/Active Ingredients and Dosage

PRODUCT	MAIN INGREDIENTS	*DOSAGE INFORMATION
Lysoform® (RKI)	Formaldehyd / anion.Tensides	1.0 %
Lysoform® Fuß-Spray	Propan-1-ol/Menthol/ Frischpflanzenextrakt/Lactol	Ready for use
Lysoformin® 3000 (VAH, IHO)	Glyoxal/glutaraldehyde/Quat/ anion. Tenside	Surfaces 0.25 % Instruments (incl. Tbc) 0.5 % incl. Mycobact. terrae 1.5 %
Lysoformin® Plus-Foam (IHO)	Polyhexamethylene biguanide / HCl/Bis-(3-aminopropyl) dodecylamine/Quat	Ready for use
Lysoformin® rapid (VAH, IHO)	Bis-(3-aminopropyl)dodecyl- amine/Quat	0.5 % – 1.0 %
Lysoformin® rapid wipes (VAH)	Bis-(3-aminopropyl)dodecyl- amine/Quat	Ready for use
Lysoformin® spezial (VAH, IHO)	Polyhexamethylene biguanide - HCl/ Quat / nichtion. Tenside	0.25 %
Majola®-H5-Creme	Oleyoleate /Fatty alcohol	Ready for use
Neoseptin® (VAH, IHO)	Ethanol	Ready for use
Trichloro® (VAH, RKI, IHO)	Tenside/caring agents	0.5 %
Wasa® 2000	Tenside / caring agents	Ready for use
Wasa®-Reinigungsschaum	anion.Tenside/refattening agents/caring agents	Ready for use
Wasa®-Soft	Tenside/refattening agents	Ready for use

RKI = Robert-Koch-Institut listing of tested and certified disinfectants
 VAH = Verbund für Angewandte Hygiene e.V. – listing of disinfectants
 IHO = Industrieverband Hygiene und Oberflächenschutz – listing of disinfectants
 DVG = Deutsche Veterinärmedizinische Gesellschaft e.V. – Listing for food and livestock industries

* Read product information carefully!

HyBag®-System



Free of perfumes and colourants

- ACO®-care sensitive**
- Skin care sensitive

Art.-No. **Shipping Unit/Container**
0738 E 10 x 700 ml-bag



Free of perfumes and colourants

- ACO®-derm sensitive**
- Alcoholic hand disinfectant

Art.-No. **Shipping Unit/Container**
0706 E 10 x 700 ml-bag

Use disinfectants safely! Always read label and product information before use.

German BAuA Reg.-No.: N-13323



Free of perfumes and colourants

- ACO®-san sensitive**
- Mild washing lotion

Art.-No. **Shipping Unit/Container**
0747 E 10 x 700 ml-bag



- AHD 2000® G**
- Hand disinfection

Art.-No. **Shipping Unit/Container**
5648 E 10 x 700 ml-bag

Use disinfectants safely! Always read label and product information before use.
German BAuA Reg.-No.: N-85424



LYSOFORM HYBAG®-DISPENSER

Set Lysoform HyBag®-dispenser + drip tray

- Drip tray
- Bracket
- Single-use pump

Art.-No.
5851
5852
5853
4743

Dispenser and bag are sold separately.

Product Overview from A to Z



ACO®-care sensitive

- Skin care sensitive

Art.-No.	Shipping Unit/Container
0736 E	12 x 500 ml-bottle
0738 E	10 x 700 ml-bag



ACO®-derm sensitive

- Alcoholic hand disinfectant

Art.-No.	Shipping Unit/Container
0701 E	20 x 125 ml-bottle
0705 E	12 x 500 ml-bottle
0706 E	10 x 700 ml-bag
0707 E	6 x 1 litre-bottle
0709 E	2 x 5 litre-canister

Use disinfectants safely! Always read label and product information before use.

German BAuA Reg.-No.: N-91332



ACO®-san sensitive

- Skin care sensitive

Art.-No.	Shipping Unit/Container
0745 E	12 x 500 ml-bottle
0747 E	10 x 700 ml-bag



AERODESIN®

- Rapid disinfection/Wiping disinfection of medical devices
(Surface disinfection acc. to VAH-methods)

Art.-No.	Shipping Unit/Container
9896 E	12 x 1 litre-bottle

Use disinfectants safely! Always read label and product information before use.

German BAuA Reg.-No.: N-12735, N-12737

CE 0482

Product Overview from A to Z



AERODESIN® 2000

CE⁰⁴⁸²

- Rapid disinfection of medical devices
(Surface disinfection acc. to VAH-methods)

Art.-No.	Shipping Unit/Container
9736 E	12 x 1 litre-bottle with spray head
9725 E	2 x 5 litre-canister

Use disinfectants safely! Always read label and product information before use.

German BAuA Reg.-No.: N-12635, N-12630



AHD 2000® G

- Hand disinfection acc. to VAH-methods

Art.-No.	Shipping Unit/Container
5576 E	20 x 125 ml-bottle
5602 E	15 x 250 ml-bottle with Spray pump
8505 E	12 x 500 ml-bottle
5584 E	6 x 1 litre-bottle
9875 E	2 x 5 litre-canister

Use biocides safely! Always read label and product information before use.
German BAuA Reg.-No.: N-109526



AHD 2000®/CH (for Sterisol Dispenser)

- Hand disinfection acc. to VAH-methods

Art.-No.	Shipping Unit/Container
5612	12 x 700 ml-bag

Use biocides safely! Always read label and product information before use.
German BAuA Reg.-No.: N-85424



AHD 2000® G (HyBag-System)

- Hand disinfection

Art.-No.	Shipping Unit/Container
5648 E	10 x 700 ml-bag

Use biocides safely! Always read label and product information before use.
German BAuA Reg.-No.: N-85424



ALDASAN® 2000 / ALDOSAN® 2000 – CH

CE⁰⁴⁸²

- Disinfection of medical devices by wiping or submersion
(Surface and instrument disinfection acc. to VAH-methods)

Art.-No.	Shipping Unit/Container
6241 E	6 x 1 litre-bottle
6255 E	2 x 5 litre-canister

Use disinfectants safely! Always read label and product information before use.

Product Overview from A to Z



ALMYROL®

CE 0482

- Disinfection of medical devices by wiping or submersion (instrument disinfection acc. to VAH-methods)

Art.-No.	Shipping Unit/Container
7214 E	6 x 1 litre-bottle
7205 E	2 x 5 litre-canister

Use disinfectants safely! Always read label and product information before use.



AMOSEPT®

Zul.-Nr.: 6339951.00.00

- Hygienic Hand wash and cleaning acc. to VAH-methods

Art.-No.	Shipping Unit/Container
4755 E	20 x 500 ml-bottle
4778 E	12 x 1 litre-bottle

Use biocides safely! Always read label and product information before use.



BLANISOL®-PUR

CE

- Cleaning of instruments and surfaces of medical devices

Art.-No.	Shipping Unit/Container
2113 E	12 x 1 litre-bottle

Use medical devices safely! Always read label and product information before use.



BLANISOL® SPEZIALSCHAUM

- Kitchen- and bath cleaner

Art.-No.	Shipping Unit/Container
0501 E	12 x 1 litre-bottle
0509 E	2 x 5 litre-canister

Use cleaning agents safely! Always read label and product information before use.

Product Overview from A to Z



CLORINA®

CE 0482

- Wiping disinfection of medical devices
(Surface disinfection acc. to VAH-methods)

Art.-No.	Shipping Unit/Container
1660 E	10 x 640 g container
1678 E	1 x 4,5 kg container

Use disinfectants safely! Always read label and product information before use.

German BAuA Reg.-No.: N-12730, N-12731, N-12732, N-12733



FUGATEN®

CE 0482

- Wiping disinfection of medical devices
(Surface disinfection acc. to VAH-methods)

Art.-No.	Shipping Unit/Container
9685 E	12 x 1 litre-bottle
2299 E	12 x 1 litre-bottle unscented
2058 E	10 x 500 ml-bottle incl. spray head
2292 E	12 x 1 litre-bottle incl. spray head
2272 E	12 x 1 litre-bottle incl. spray head unscented
9825 E	2 x 5 litre-canister
9823 E	2 x 5 litre-canister unscented

Use disinfectants safely! Always read label and product information before use.

German BAuA Reg.-No.: N-14992, N-14994



HOSPISEPT®

Zul.-Nr.: 6339833.00.00

- Hand and Skin disinfection by mechanical action acc. to VAH

Art.-No.	Shipping Unit/Container
9416 E	20 x 125 ml-bottle
9654 E	15 x 250 ml-bottle with spray head
2571 E	12 x 500 ml-bottle
2567 E	12 x 1 litre bottle
2605 E	2 x 5 litre canister

Use disinfectants safely! Always read label and product information before use.



HOSPISEPT®-WIPES

CE 0482

- Wiping disinfection of medical devices

Art.-No.	Shipping Unit/Container
3595 E	1 Box containing 15 single wipes
3585 E	6 containers (100 wipes each)
3575 E	6 refill bags (100 wipes each)

Use disinfectants safely! Always read label and product information before use.

German BAuA Reg.-No.: N-12666, N-12667, N-12670

Product Overview from A to Z



LUPHENIL®

- Intensive skin care and skin regeneration

Art.-No. **Shipping Unit/Container**

3301 E 20 x 50 ml-tube



LY®-DRY

- Rapid binder for liquids

Art.-No. **Shipping Unit/Container**

8871 E 12 x 400 g-container



LY®-FRESH

- odour eliminator

Art.-Nr. **Shipping Unit/Container**

8872 E 20 x 200 ml-bottle incl. spray head



LYISOFORM®

- Disinfection of medical devices by wiping and of laundry

Art.-No. **Shipping Unit/Container**

6545 E 20 x 500 ml-bottle

6566 E 12 x 1 litre-bottle

6444 E 2 x 5 litre-canister

CE 0482

Use disinfectants safely! Always read label and product information before use.

Product Overview from A to Z



LYSOFORM® FUSS-SPRAY

Art.-No.	Shipping Unit/Container
1012 E	15 x 250 ml-bottle with spray pump



LYSOFORMIN® 3000

CE 0482

- Disinfection of medical devices by wiping or submersion (Disinfection of surfaces and instruments acc. to VAH-methods)

Art.-No.	Shipping Unit/Container
9508 E	200 x 20 ml-dosage-bag = 4 litre per shipping unit
9256 E	12 x 1 litre-bottle
9265 E	2 x 5 litre-canister
9330 E	1 x 20 litre-canister

Use disinfectants safely! Always read label and product information before use.
German BAuA Reg.-No.: N-12657, N-12658, N-12659



LYSOFORMIN® PLUS-SCHAUM

CE 0482

- Spray disinfection of medical devices (surface disinfection)

Art.-No.	Shipping Unit/Container
2505 E	12 x 1 litre-bottle with foam head
2509 E	2 x 5 litre-canister

Use disinfectants safely! Always read label and product information before use.

German BAuA Reg.-No.: N-76300



LYSOFORMIN® RAPID

CE 0482

- Wiping disinfection of medical devices (Surface disinfection acc. to VAH-methods)

Art.-No.	Shipping Unit/Container
9238 E	12 x 1 litre-bottle
9240 E	2 x 5 litre-canister

Use disinfectants safely! Always read label and product information before use.

German BAuA Reg.-No.: N-70957

Product Overview from A to Z



LYSOFORMIN® RAPID WIPES

CE 0482

- Impregnated wipes for the wiping disinfection of medical devices (Surface disinfection acc. to VAH-methods)

Art.-No.	Shipping Unit/Container
9235 E	4 Bags with 100-wipes rolls 1 Extraction cap 1 Mini cutting knife

Use Lysoformin rapid wipes (contains disinfectant) safely! Use disinfectants safely! Always read label and product information before use.

German BAuA Reg.-No.: N-73365



LYSOFORMIN® SPEZIAL

CE 0482

- Wiping disinfection of medical devices (Surface disinfection acc. to VAH-methods)

Art.-No.	Shipping Unit/Container
9287 E	200 x 40 ml-Dosage-bag = 8 litre per shipping unit
9291 E	12 x 1 litre-bottle
9325 E	2 x 5 litre-canister
9320 E	1 x 20 litre-canister

Use disinfectants safely! Always read label and product information before use.

German BAuA Reg.-No.: N-76302



MAJOLA® H5-CREME

- Skin care sensitive

Art.-No.	Shipping Unit/Container
2312 E	10 x 100 ml-tube
9794 E	12 x 500 ml-bottle
9841 E	6 x 1 litre-bottle



NEOSEPTIN®

CE 0482

- Hand disinfection acc. to VAH-methods and disinfection of medical devices

Art.-No.	Shipping Unit/Container
2700 E	12 x 500 ml-bottle
2705 E	6 x 1 litre-bottle

Use disinfectants safely! Always read label and product information before use.

German BAuA Reg.-No.: N-39624, N-39623, N-39625

Product Overview from A to Z



TRICHLOROL®

CE 0482

- Wiping disinfection of medical devices
(Surface disinfection acc. to VAH-methods)

Art.-No. Shipping Unit/Container

1701 E	10 x 500 g-container
1720 E	1 x 5 kg-container

Use disinfectants safely! Always read label and product information before use.

German BAuA Reg.-No.: N-12753, N-12754, N-12755



WASA® 2000

- Washing and care lotion

Art.-No. Shipping Unit/Container

6104 E	20 x 125 ml-bottle
8498 E	20 x 500 ml-bottle
6118 E	12 x 1 litre-bottle
5545 E	2 x 5 litre-canister



WASA®-REINIGUNGSSCHAUM

- Free of greenhouse gas emissions

Art.-No. Shipping Unit/Container

1835 E	18 x 200 ml-bottle with foamhead
4825 E	2 x 5 liter-canister



WASA®-SOFT

- Washing and caring lotion

Art.-No. Shipping Unit/Container

4845 E	20 x 500 ml-bottle
4875 E	12 x 1 litre-bottle
4895 E	2 x 5 litre-canister

Accessories



LYSOFORM®-SOFT-VLIES

Art.-No.

- Universal fleece cloth dispenser for use with aqueous Disinfectants or cleaning agents

Shipping Unit/Container

6 dispenser-containers 3400
6 dry rolls (100 wipes each) 3402

EMPTY BOTTLE (1 litre)

Art.-No.

30110

ORIGINALITÄTS-KSV CAP 24 NATUR 2 MM FOR EMPTY BOTTLES

30420

EMPTY SPRAY BOTTLE (1 litre)

30598

SPRAY HEAD FOR 1-L-SPRAY BOTTLE

30132




DOSAGE HEAD for 1000 ml bottle

Art.-No.
9774

ADAPTER

9776



DOSAGE AID DH 1

Art.-No.
9170

DOSAGE AID DH 91

9174

Please Note when ordering:

1. Dosage Aid 91 includes all available accessories.
2. The product with which the dosage aid will be applied should be noted in the order.



DOSAGE PUMP for 5-litre canisters

Art.-No.

10 ml per pump
20 ml per pump

3134
3135



MEASURING CUP 100 ml

Art.-No.
4600

Accessories

	2-NECK-DOSING BOTTLE 1 litre	Art.-No. 9773
	DISPENSER for 500 ml bottle for 1 litre-bottle	Art.-No. 9771 9772
	LYSOFORM HYBAG®-DISPENSER Set Lysoform HyBag®-DISPENSER + drip tray Drip tray Bracket Single-use pump	Art.-No. 5851 5852 5853 4743
	Dispenser and bag sold separately.	
	LYDOS F WALL-MOUNTED DISPENSER for 500 ml bottle LYDOS 1000 WALL-MOUNTED DISPENSER for 1 litre-bottle	Art.-No. 5810 5763
	CANISTER KEY for 5 litre-canister	Art.-No. 7073
	TRANSFER TAP for 5 litre-canister	Art.-No. 3156
	LYSOFORM DISPENSER COLUMN Brochure holder for Lysoform Dispenser Column	Art.-No. 2890 2891

Dosing Table

Lysoform Dr. Hans Rosemann GmbH
www.lysoform.com



amount of utility dilution \ utility concentration	0.25%	0.50%	0.75%	1.00%	1.50%	2.00%	3.00%	4.00%	5.00%
1 liter	2.5 ml	5 ml	7.5 ml	10 ml	15 ml	20 ml	30 ml	40 ml	50 ml
2 liter	5.0 ml	10 ml	15.0 ml	20 ml	30 ml	40 ml	60 ml	80 ml	100 ml
3 liter	7.5 ml	15 ml	22.5 ml	30 ml	45 ml	60 ml	90 ml	120 ml	150 ml
4 liter	10.0 ml	20 ml	30.0 ml	40 ml	60 ml	80 ml	120 ml	160 ml	200 ml
5 liter	12.5 ml	25 ml	37.5 ml	50 ml	75 ml	100 ml	150 ml	200 ml	250 ml
6 liter	15.0 ml	30 ml	45.0 ml	60 ml	90 ml	120 ml	180 ml	240 ml	300 ml
7 liter	17.5 ml	35 ml	52.5 ml	70 ml	105 ml	140 ml	210 ml	280 ml	350 ml
8 liter	20.0 ml	40 ml	60.0 ml	80 ml	120 ml	160 ml	240 ml	320 ml	400 ml
9 liter	22.5 ml	45 ml	67.5 ml	90 ml	135 ml	180 ml	270 ml	360 ml	450 ml
10 liter	25.0 ml	50 ml	75.0 ml	100 ml	150 ml	200 ml	300 ml	400 ml	500 ml
15 liter	37.5 ml	75 ml	112.5 ml	150 ml	225 ml	300 ml	450 ml	600 ml	750 ml
20 liter	50.0 ml	100 ml	150.5 ml	200 ml	300 ml	400 ml	600 ml	800 ml	1.0 liter
25 liter	62.5 ml	125 ml	187.5 ml	250 ml	375 ml	500 ml	750 ml	1.0 liter	1.25 liter
30 liter	75.0 ml	150 ml	225.0 ml	300 ml	450 ml	600 ml	900 ml	1.2 liter	1.5 liter
40 liter	100.0 ml	200 ml	300.0 ml	400 ml	600 ml	800 ml	1.2 liter	1.6 liter	2.0 liter
50 liter	125.0 ml	250 ml	375.0 ml	500 ml	750 ml	1.0 liter	1.5 liter	2.0 liter	2.5 liter

Important: Always fill with water (approx. room temperature) first and then add the disinfectant!

The amount of water is calculated by subtracting the amount of disinfectant from the required amount of the use dilution.






Gaps in coverage after hand disinfection

The areas marked in orange represents common gaps in coverage after hand disinfection



The six steps of hand disinfection:

Standard hygienic hand rub: A sufficient amount of hand disinfectant should be applied to the cupped, dry palm, and be rubbed in according to the instructions below, for a minimum of 30 seconds. Disinfectant should be spread without gaps down to the wrists to ensure thorough moistening of the hands. The movements described in each of the six steps should be performed five times, before moving on to the next step. Caution: Ensure thorough moistening of the hands throughout the entire application period.





LYSOFORM
Seit ▪ Since ▪ Depuis 1900

Lysoform Dr. Hans Rosemann GmbH
Kaiser-Wilhelm-Str. 133, D-12247 Berlin
Telefon +49 30 77992-0, Telefax +49 30 77992-269
www.lysoform.com, e-mail: kontakt@lysoform.de